



# SLC10A6 Antibody, FITC conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA665974LC01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                 |
| <b>Uniprot No.</b>         | Q3KNW5  |
| <b>Immunogen</b>           | Recombinant Human Solute carrier family 10 member 6 protein (312-377AA)                       |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | FITC  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4              |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Solute carrier family 10 member 6 (Sodium-dependent organic anion transporter), SLC10A6, SOAT |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Research Area</b>       | Others  |
| <b>Target Names</b>        | SLC10A6   |