



# KRTAP13-4 Antibody, Biotin conjugated

Product Code	CSB-PA667309LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3LI77
Immunogen	Recombinant Human Keratin-associated protein 13-4 protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin-associated protein 13-4, KRTAP13-4, KAP13.4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KRTAP13-4