



TTLL9 Antibody, Biotin conjugated

Product Code	CSB-PA667484LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3SXZ7
Immunogen	Recombinant Human Probable tubulin polyglutamylase TTLL9 protein (1-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable tubulin polyglutamylase TTLL9 (EC 6.-.-.-) (Tubulin--tyrosine ligase-like protein 9), TTLL9, C20orf125
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TTLL9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.