



SP6 Antibody, Biotin conjugated

Product Code	CSB-PA668784LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3SY56
Immunogen	Recombinant Human Transcription factor Sp6 protein (101-376AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor Sp6 (Krueppel-like factor 14), SP6, KLF14
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.