



SLC36A2 Antibody, Biotin conjugated

Product Code	CSB-PA672077LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q495M3
Immunogen	Recombinant Human Proton-coupled amino acid transporter 2 protein (1-58AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proton-coupled amino acid transporter 2 (Proton/amino acid transporter 2) (Solute carrier family 36 member 2) (Tramdorin-1), SLC36A2, PAT2 TRAMD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC36A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.