



TIGIT Antibody, Biotin conjugated

Product Code	CSB-PA675446LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q495A1
Immunogen	Recombinant Human T-cell immunoreceptor with Ig and ITIM domains protein (22-141AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell immunoreceptor with Ig and ITIM domains (V-set and immunoglobulin domain-containing protein 9) (V-set and transmembrane domain-containing protein 3), TIGIT, VSIG9 VSTM3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TIGIT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.