



USH1G Antibody, Biotin conjugated

Product Code	CSB-PA676572LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q495M9
Immunogen	Recombinant Human Usher syndrome type-1G protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Usher syndrome type-1G protein (Scaffold protein containing ankyrin repeats and SAM domain), USH1G, SANS
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	USH1G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.