



USH1G Antibody, Biotin conjugated

Product Code	CSB-PA676572LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q495M9
Immunogen	Recombinant Human Usher syndrome type-1G protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Usher syndrome type-1G protein (Scaffold protein containing ankyrin repeats and SAM domain), USH1G, SANS
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	USH1G