




GTDC1 Antibody

Product Code	CSB-PA676737LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4AE62
Immunogen	Recombinant Human Glycosyltransferase-like domain-containing protein 1 protein (1-292AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycosyltransferase-like domain-containing protein 1 (Mannosyltransferase-like protein Xa) (Mat-Xa), GTDC1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GTDC1

Image



Western blot

All lanes: GTDC1 antibody at 4µg/ml + A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 53, 44, 34, 50, 38 kDa

Observed band size: 53 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.