



ssaA1 Antibody, FITC conjugated

| | |
|----------------------------|---------------------------------------------------------------------------------------------|
| Product Code | CSB-PA683042LC01FLB |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q5HCY4 |
| Immunogen | Recombinant Staphylococcus aureus Staphylococcal secretory antigen ssaA1 protein (27-255AA) |
| Raised In | Rabbit |
| Species Reactivity | Staphylococcus aureus |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Staphylococcal secretory antigen ssaA1, ssaA1 |
| Immunogen Species | Staphylococcus aureus (strain COL) |
| Research Area | Others |
| Target Names | ssaA1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |