



ssaA1 Antibody, Biotin conjugated

Product Code	CSB-PA683042LD01FLB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5HCY4
Immunogen	Recombinant Staphylococcus aureus Staphylococcal secretory antigen ssaA1 protein (27-255AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Staphylococcal secretory antigen ssaA1, ssaA1
Immunogen Species	Staphylococcus aureus (strain COL)
Research Area	Others
Target Names	ssaA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.