





TTC32 Antibody

Product Code	CSB-PA686303LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5I0X7
Immunogen	Recombinant Human Tetratricopeptide repeat protein 32 protein (1-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 32 (TPR repeat protein 32), TTC32
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TTC32

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA686303LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA686303LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.