



TTC32 Antibody, HRP conjugated

Product Code	CSB-PA686303LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5I0X7
Immunogen	Recombinant Human Tetratricopeptide repeat protein 32 protein (1-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 32 (TPR repeat protein 32), TTC32
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TTC32
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.