



MRPL51 Antibody, HRP conjugated

Product Code	CSB-PA687045LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4U2R6
Immunogen	Recombinant Human 39S ribosomal protein L51, mitochondrial protein (32-128AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L51, mitochondrial (L51mt) (MRP-L51) (Mitochondrial large ribosomal subunit protein mL51) (bMRP-64) (bMRP64), MRPL51, MRP64
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRPL51
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.