



Ribosome-inactivating protein gynostemmin Antibody, Biotin conjugated

Product Code	CSB-PA688696LD01GAAZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5F0I3
Immunogen	Recombinant Gynostemma pentaphyllum Ribosome-inactivating protein gynostemmin protein (24-300AA)
Raised In	Rabbit
Species Reactivity	Gynostemma pentaphyllum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosome-inactivating protein gynostemmin (EC 3.2.2.22) (rRNA N-glycosidase)
Immunogen Species	Gynostemma pentaphyllum (Jiaogulan) (Sweet tea vine)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.