



TMEM67 Antibody, Biotin conjugated

Product Code	CSB-PA689243LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5HYA8
Immunogen	Recombinant Human Meckelin protein (241-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Meckelin (Meckel syndrome type 3 protein) (Transmembrane protein 67), TMEM67, MKS3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TMEM67
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.