



GCSAML Antibody, Biotin conjugated

Product Code	CSB-PA689332LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JQS6
Immunogen	Recombinant Human Germinal center-associated signaling and motility-like protein (1-135AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Germinal center-associated signaling and motility-like protein, GCSAML, C1orf150
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GCSAML
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.