



fabF Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA692016LB01FLB |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q5HHA1 |
| Immunogen | Recombinant Staphylococcus aureus 3-oxoacyl-[acyl-carrier-protein] synthase 2 protein (1-414AA) |
| Raised In | Rabbit |
| Species Reactivity | Staphylococcus aureus |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | 3-oxoacyl-[acyl-carrier-protein] synthase 2 (EC 2.3.1.179) (3-oxoacyl-[acyl-carrier-protein] synthase II) (Beta-ketoacyl-ACP synthase II) (KAS II), fabF |
| Immunogen Species | Staphylococcus aureus (strain COL) |
| Research Area | Others |
| Target Names | fabF |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |