



GLE1 Antibody, FITC conjugated

Product Code	CSB-PA700651EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53GS7
Immunogen	Recombinant Human Nucleoporin GLE1 protein (141-380AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoporin GLE1 (hGLE1) (GLE1-like protein), GLE1, GLE1L
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GLE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.