





BEX5 Antibody

Product Code	CSB-PA705160LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5H9J7
Immunogen	Recombinant Human Protein BEX5 protein (1-111AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein BEX5 (Brain-expressed X-linked protein 5) (NGFRAP1-like protein 1) (Nerve growth factor receptor-associated protein 2), BEX5, NADE2, NGFRAP1L1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	BEX5

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA705160LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA705160LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.