



MAP3K19 Antibody, HRP conjugated

Product Code	CSB-PA707114LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q56UN5
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 19 protein (45-168AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 19 (EC 2.7.11.1) (Regulated in COPD, protein kinase) (SPS1/STE20-related protein kinase YSK4), MAP3K19, RCK YSK4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAP3K19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.