



MAP3K19 Antibody, FITC conjugated

Product Code	CSB-PA707114LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q56UN5
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 19 protein (45-168AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 19 (EC 2.7.11.1) (Regulated in COPD, protein kinase) (SPS1/STE20-related protein kinase YSK4), MAP3K19, RCK YSK4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAP3K19