



PAGE2B Antibody, FITC conjugated

Product Code	CSB-PA708481LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JRK9
Immunogen	Recombinant Human Putative G antigen family E member 3 protein (1-111AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative G antigen family E member 3 (Prostate-associated gene 2B protein) (PAGE-2B), PAGE2B, GAGEE3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PAGE2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.