



TMEM255A Antibody, Biotin conjugated

Product Code	CSB-PA708482LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JRV8
Immunogen	Recombinant Human Transmembrane protein 255A protein (110-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 255A (Protein FAM70A), TMEM255A, FAM70A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TMEM255A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.