




SPIN3 Antibody

Product Code	CSB-PA708489LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JUX0
Immunogen	Recombinant Human Spindlin-3 protein (1-77AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spindlin-3 (Spindlin-like protein 3) (SPIN-3), SPIN3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPIN3

Image



Western blot

All lanes: SPIN3 antibody at 6µg/ml

Lane 1: A549 whole cell lysate

Lane 2: HeLa whole cell lysate

Lane 3: A431 whole cell lysate

Lane 4: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 30, 9 kDa

Observed band size: 30, 70 kDa