



SPIN3 Antibody, Biotin conjugated

Product Code	CSB-PA708489LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JUX0
Immunogen	Recombinant Human Spindlin-3 protein (1-77AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spindlin-3 (Spindlin-like protein 3) (SPIN-3), SPIN3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPIN3