



KHDC1 Antibody, Biotin conjugated

Product Code	CSB-PA709298LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4VXA5
Immunogen	Recombinant Human KH homology domain-containing protein 1 protein (74-237AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	KH homology domain-containing protein 1, KHDC1, C6orf147 C6orf148
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	KHDC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.