



SH3D19 Antibody, FITC conjugated

Product Code	CSB-PA711356LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5HYK7
Immunogen	Recombinant Human SH3 domain-containing protein 19 protein (468-767AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 domain-containing protein 19 (ADAM-binding protein Eve-1) (EEN-binding protein) (EBP), SH3D19
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SH3D19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.