



MIA3 Antibody

Product Code	CSB-PA711448LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JRA6
Immunogen	Recombinant Human Transport and Golgi organization protein 1 homolog protein (1447-1626AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transport and Golgi organization protein 1 homolog (TANGO1) (C219-reactive peptide) (D320) (Melanoma inhibitory activity protein 3), MIA3, KIAA0268 TANGO
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MIA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.