



MIA3 Antibody, Biotin conjugated

Product Code	CSB-PA711448LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JRA6
Immunogen	Recombinant Human Transport and Golgi organization protein 1 homolog protein (1447-1626AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transport and Golgi organization protein 1 homolog (TANGO1) (C219-reactive peptide) (D320) (Melanoma inhibitory activity protein 3), MIA3, KIAA0268 TANGO
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MIA3