



FAM78A Antibody, HRP conjugated

Product Code	CSB-PA711454LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JUQ0
Immunogen	Recombinant Human Protein FAM78A protein (21-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM78A, FAM78A, C9orf59
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM78A