



FAM78A Antibody, FITC conjugated

Product Code	CSB-PA711454LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JUQ0
Immunogen	Recombinant Human Protein FAM78A protein (21-283AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM78A, FAM78A, C9orf59
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM78A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.