



KRBOX4 Antibody, Biotin conjugated

Product Code	CSB-PA711455LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JUW0
Immunogen	Recombinant Human KRAB domain-containing protein 4 protein (1-98AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	KRAB domain-containing protein 4 (KRAB box domain-containing protein 4), KRBOX4, ZNF673
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	KRBOX4