



KRBOX4 Antibody, Biotin conjugated

Product Code	CSB-PA711455LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JUW0
Immunogen	Recombinant Human KRAB domain-containing protein 4 protein (1-98AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	KRAB domain-containing protein 4 (KRAB box domain-containing protein 4), KRBOX4, ZNF673
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	KRBOX4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.