



C6orf136 Antibody, HRP conjugated

Product Code	CSB-PA713102LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5SQH8
Immunogen	Recombinant Human Uncharacterized protein C6orf136 protein (1-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C6orf136, C6orf136
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C6orf136
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.