




C9orf64 Antibody

Product Code	CSB-PA713128LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T6V5
Immunogen	Recombinant Human Queuosine salvage protein (182-341AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Queuosine salvage protein, C9orf64
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C9orf64

Image



Western blot

All lanes: C9orf64 antibody at 1μg/ml

Lane 1: Hela whole cell lysate

Lane 2: K562 whole cell lysate


Lane 3: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 40 kDa


Observed band size: 40 kDa



Immunohistochemistry of paraffin-embedded

human liver cancer using CSB-

PA713128LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA713128LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.