



GLT8D1 Antibody, FITC conjugated

Product Code	CSB-PA715010LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68CQ7
Immunogen	Recombinant Human Glycosyltransferase 8 domain-containing protein 1 protein (33-332AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycosyltransferase 8 domain-containing protein 1 (EC 2.4.1.-), GLT8D1, GALA4A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GLT8D1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.