



KIF24 Antibody, Biotin conjugated

Product Code	CSB-PA716354LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T7B8
Immunogen	Recombinant Human Kinesin-like protein KIF24 protein (1168-1368AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF24, KIF24, C9orf48
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KIF24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.