





KRT40 Antibody

Product Code	CSB-PA718197LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6A162
Immunogen	Recombinant Human Keratin, type I cytoskeletal 40 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cytoskeletal 40 (Cytokeratin-40) (CK-40) (Keratin-40) (K40) (Type I hair keratin Ka36), KRT40, KA36
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT40

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA718197LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA718197LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.