



KRT40 Antibody, HRP conjugated

Product Code	CSB-PA718197LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6A162
Immunogen	Recombinant Human Keratin, type I cytoskeletal 40 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cytoskeletal 40 (Cytokeratin-40) (CK-40) (Keratin-40) (K40) (Type I hair keratin Ka36), KRT40, KA36
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT40
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.