



FCRL6 Antibody, Biotin conjugated

Product Code	CSB-PA718779LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6DN72
Immunogen	Recombinant Human Fc receptor-like protein 6 protein (20-307AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fc receptor-like protein 6 (FcR-like protein 6) (FcRL6) (Fc receptor homolog 6) (FcRH6) (IFGP6), FCRL6, FCRH6
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCRL6