



FCRL6 Antibody, Biotin conjugated

Product Code	CSB-PA718779LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6DN72
Immunogen	Recombinant Human Fc receptor-like protein 6 protein (20-307AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fc receptor-like protein 6 (FcR-like protein 6) (FcRL6) (Fc receptor homolog 6) (FcRH6) (IFGP6), FCRL6, FCRH6
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCRL6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.