




LCE2A Antibody

Product Code	CSB-PA719417LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TA79
Immunogen	Recombinant Human Late cornified envelope protein 2A protein (1-106AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Late cornified envelope protein 2A (Late envelope protein 9), LCE2A, LEP9
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LCE2A

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA719417LA01HU at dilution of 1:100