




KRT79 Antibody

Product Code	CSB-PA719894LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5XKE5
Immunogen	Recombinant Human Keratin, type II cytoskeletal 79 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type II cytoskeletal 79 (Cytokeratin-79) (CK-79) (Keratin-6-like) (Keratin-6L) (Keratin-79) (K79) (Type-II keratin Kb38), KRT79, K6L KB38 KRT6L
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	KRT79

Image



Western blot

All lanes: KRT79 antibody at 1µg/ml

Lane 1: Mouse lung tissue

Lane 2: A549 whole cell lysate

Lane 3: 293T whole cell lysate

Lane 4: Mouse gonadal tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 58 kDa

Observed band size: 58 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.