



# KRT79 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA719894LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q5XKE5
<b>Immunogen</b>	Recombinant Human Keratin, type II cytoskeletal 79 protein (1-300AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Keratin, type II cytoskeletal 79 (Cytokeratin-79) (CK-79) (Keratin-6-like) (Keratin-6L) (Keratin-79) (K79) (Type-II keratin Kb38), KRT79, K6L KB38 KRT6L
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Tags & Cell Markers
<b>Target Names</b>	KRT79
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.