



SCARA3 Antibody, HRP conjugated

Product Code	CSB-PA721489LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6AZY7
Immunogen	Recombinant Human Scavenger receptor class A member 3 protein (201-466AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Scavenger receptor class A member 3 (Cellular stress response gene protein), SCARA3, CSR
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SCARA3