



SCARA3 Antibody, Biotin conjugated

Product Code	CSB-PA721489LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6AZY7
Immunogen	Recombinant Human Scavenger receptor class A member 3 protein (201-466AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Scavenger receptor class A member 3 (Cellular stress response gene protein), SCARA3, CSR
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SCARA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.