



PVRIG Antibody, Biotin conjugated

Product Code	CSB-PA721905LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6DKI7
Immunogen	Recombinant Human Transmembrane protein PVRIG protein (193-326AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein PVRIG (CD112 receptor) (CD112R) (Poliovirus receptor-related immunoglobulin domain-containing protein), PVRIG, C7orf15
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PVRIG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.