



ZDHHC20 Antibody, Biotin conjugated

Product Code	CSB-PA722800LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5W0Z9
Immunogen	Recombinant Human Probable palmitoyltransferase ZDHHC20 protein (75-130AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable palmitoyltransferase ZDHHC20 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 20) (DHHC-20), ZDHHC20
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZDHHC20
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.