




EGFLAM Antibody

Product Code	CSB-PA723649LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q63HQ2
Immunogen	Recombinant Human Pikachurin protein (25-300AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pikachurin (Agrin-like protein) (EGF-like, fibronectin type-III and laminin G-like domain-containing protein), EGFLAM, AGRINL AGRNL PIKA
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	EGFLAM

Image



Western blot
All lanes: EGFLAM antibody at 3µg/ml + Mouse heart tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 112, 111, 42, 85, 17 kDa
Observed band size: 112 kDa